IDENCIFICACION OF A RAPID STD TESC

Kehinde Sogbesan

Breaking News "Std's have been on the rise for Sixth consecutive years in the United States"

CDC Reports as of 2021


STD EPIDEMIC AC LARGE

STD- sexually transmitted disease

- 1 million STDs acquired every day
- Chlamydia trachomatis and Neisseria gonorrhoea topping the list annually
- 1 in 5 people have an STD

Symptoms and Consequences of Chlamydia and Gonorrhea:

- Painful Urination
- Discharge in men and women
- Infertility
- Infections to Newborn
- Increased risk of HIV/AIDS
- Pelvic Inflammatory Disease



Chlamy dia in Uterus

RAPID POINC OF Care Tescing

RAPID POINT OF CARE TESTING

- POC- point of care diagnostic tests one can administer on themselves
- Lateral Flow Devices are common POC Tests

C Abbott

Examples:

- Pregnancy Test
- At Home COVID-Test
- HIV and Syphilis Finger Poke Test

Faster Detection \rightarrow Faster Intervention \rightarrow Decrease in STD Rates

1^{SE} RAPID POC CHLAMYDIA and Gonorrhea test

<u>Sciencific</u>

- What part of the bacterial Infections to target?
- What is the procedure to target them?
- How can we reproduce this study in urine?

<u>Encrepreneurial</u>

- Who are my competitors in the industry?
- What are the current challenges in the industry?
- Who are the stakeholders and what is the best implementation plan for the product?

Encrepreneurial Lens





ICare Fluid LFD

Most Competitive

- Urethra Swab and Endocervix
 Swab LFD
- Results in 10-15 Minutes

FIErce compecicion



BINX HEALCH ASSAY

Competitive

- Rapid PCR Test
- Results in 30 Minutes



EVERLYWELL HOME TEST

Least Competitive

- Urine samples sent to labs
- Results in 5-10 days

EVERLYWELL HOME TEST RAPHD



BINX HEALCH ASSAY AFFORDABLE AC-HOME

O)



Detects Chlamydia (Chlamydia Trachomatis) Antigen High Clinical Accuracy Reliable & Fast Results in <15 mins Tested by Certified Laboratories Worldwide ISO 13485 & GMP Certifications Easy to Use

Centified

Product of JAL Innovation Singapore

Store in a Cool Place Between 4⁰ To 30⁰ C Do Not Place Under Direct Sunlight

ICARE FLUID LFD EASY CO USE

Value proposicions



SCAKeHOLDers

A Product for the people

Identified Stakeholders:

- Customers: college students, disenfranchised communities, minorities, and sexually active individuals
- Indirect/External Stakeholders: Health Clinics, College Health Departments, Governmental Health Agencies

Is this Product Valuable?

- City Wide Survey with Colleges and Health Clinics
- Surveyed Institutions:
 - Florida State University
 - Florida Agriculture and Mechanical University
 - Neighborhood Medical Center
 - Planned Parenthood
 - Florida Department of Health



Response Summary

WHY DO INDIVIDUALS refrain from ROUCINELY STD TESC?

- Lack of Sexual Health Knowledge
- Costly: sliding scale services, insurance payment, or cash payment
- Fear, Shame, and Embarrassment

"Some students feel they can never catch a STD even though they are sexually active" COULD CHIS SCD CESC BE OFFERED IN YOUR INSCICUION?

- Overwhelming Support
- Fast accurate results helps with contact tracing

"The faster we have results the faster we can help provide proper treatment"

FINANCIAL AND OUTREACH

- Most beneficial at an affordable price or free
- Offered in universities, health clinics, outreach programs, and health stores
 - Examples: Mobile Health Units, College Health Departments, Pharmacy Stores


Maximum Outreach: Approved by FDA and administered by Department of Health

IMPLEMENCACION PLAN



SCIENCIFIC LENS


IDENCIFYING OUT TATGEC



SECUTING CHAC SPECIAL BOND WICH MOMP

Initial Idea:

- Oligosaccharides-carbohydrates containing 3-10 sugar units
- Adhesins on MOMP protein will bind to oligosaccharides to infect epithelial cells
- Drawbacks: One can not prove binding will occur in a laboratory experiment




ANCIBODY AND MOMP TOGECHER FOREVER

- Antibody- proteins that bind to antigens and provide protection for immune system
- Strong affinity to their targets
- Can be expensive and binding may not occur

Purchases for each STD:


- 2 Monoclonal Antibodies: high specificity to epitope
- 1 Polyclonal Antibodies: affinity for antigen but different epitopes




Yale COVID Antibodies

ELISA Procedure


ELISA SANDWICH ASSAY



ANCIBODY COACED IN WELL




MOMP PROCEIN BINDS CO ANCIBODY



CONJUGACE ANCIBODY BINDS CO MOMP




SUBSCRACE CURNS WELL BLUE INDICACING BINDING




SCOP SOLUCION TURNS WELL YELLOW SCOPPING CHE REACCION

CHLAMYDIA




CHLAMYDIA




BeSC COMBINACION Primary Monoclonal Ancibody & Secondary Monoclonal Ancibody

GONORLHea





GONORLHea



Besc combination Primary Polyclonal Antibody & Secondary Monoclonal Antibody

WHACCOMES NEXC?





Sepmag ELISA Biosensors LFD

Kehinde's Jiggyjuice

An Extraction Buffer that will isolate the MOMP protein in the urine and increase the affinity of the protein to the antibody for LFD Applications

– IN CURRENT DEVELOPMENT –)–

ACKNOWLEDGMENCS

Jeremiah 29: 11 "For I know the plans I have for you," declares the Lord

Angela C. Richardson Award

- Supervising Advisor: Dr. Jeff Whalen
- FSU Faculty: Dr. Steven Lenhert, Dr. Kathryn Jones, Dr. Hengli Tang
- Jim Moran College of Entrepreneurship
- FSU CRE: Dave Montez

Special Thanks to Biofront Technologies: Dr. Jason Robotham, Dr. Henry Grise, Ms. Brittany Brewers

• Friends and Family





FLORIDA STATE UNIVERSITY IM MORAN COLLEGE OF ENTREPRENEURSHIP









Many Thanks!



ALTERNATIVE RESOURCES

